


Effectiveness –
combined with
stability and security

conceptnatal 



Phototherapy wrap – BiliTex



Effectiveness – combined with security and stability with our Bilitex positioning aid specifically developed for phototherapy. Particularly during therapeutic phases the little patients are prone to stretch patterns and nervous movements due to lack of muscle tone and the unbounded surroundings, which cost a lot of energy. With the Bilitex positioning aid the preemie experiences a secure and flexible environment during the therapy, which allows the infant to adopt a relaxed position. The innovative textile material features a special mesh structure with 85% permeability to the high-intensity wavelength spectrum.


In addition, the knit fabric is very soft and breathable – the desired size can be individually adjusted using a flexible closure (Velcro®). We have developed this product in collaboration with the University Hospital Tübingen. The work of the neonatal ward is oriented to the development-promoting and individualised training programmes of the NIDCAP whose pioneering guidelines have been adopted.



- + Innovative knit fabric with 85% permeability for intensive phototherapy in the wavelength range of 450-470 nm
- + Special fit with closed, elastic foot-piece ensures flexible limitation
- + Arm and shoulder section with individually adjustable loop supports orientation towards the body's centre line
- + Ease of use due to flexible Velcro® closure
- + Breathable and thermostatic
- + 5 sizes: S < 800g, M 800g - 1,200g, L 1,200g - 2,000g, XL 2,000g - 3,000g, XXL > 3,000g
- + additionally available as flexible cloth size 60cm x 60cm
- + Individually packaged / 100% PES / 10 items per pack for each size
- + Tested and certified according to Oeko-Tex®Standard 100
- + Developed and made in Germany



BiliTex cloth 60x60cm REF 18-06-100



BiliTex S	REF 05-01-100
BiliTex M	REF 05-02-100
BiliTex L	REF 05-03-100
BiliTex XL	REF 05-04-100
BiliTex XXL	REF 05-05-100
BiliTex cloth 60x60cm	REF 18-06-100

conceptnatal 

Im Michelreis 18
72124 Pliezhausen
Germany
tel +49(0)7127 9802601
fax +49(0)7127 9802602
info@conceptnatal.de
www.conceptnatal.de

